

---

## Software **FIAStudio**

### Functions

- Control of the FIA system
- Data acquisition of the FIA channels
- Presentation and evaluation of the measurement data
- Creation of result protocols of measurement series
- Data export and data archiving

The FIA system by **MLE** is controlled by the Windows-based **FIAStudio** software developed by **Dr. Herbert Steiner Consulting**.

### Control of the FIA system

The FIA system can consist of one or more channels. The control and administration of the FIA system is organized by the **FIAStudio** software. This software enables the system to be set up and analytical methods complete with start-up and shutdown routines to be created.

**FIAStudio** also manages and controls the autosampler with dilutor and organises the occupancy of the sample tray with standards, control samples and samples.

The software modules

- **FIAStudio sample preparation** for sample pre-dilution and standard preparation
- **FIAStudio PeakInspector** as a quality assurance module for checking the measurements
- **FIAStudio rights management** to manage and ensure user rights

can be optionally integrated.

The FIA system is designed for routine operation, but it is also suitable for method development.

### Data acquisition from **FIA Nexus**

- data acquisition and management for all channels
- data storage in database
- presentation of the analytical peak of the last measurement
- status indication for the connected modules
- connection between FIA system and PC via serial interface RS 232 or USB

## Presentation and evaluation of measurement data

- display of measurement data based on the selection of series, measurement date or time
- presentation of the sample peaks including zoom function
- marking of incorrect measurements as outliers
- presentation of the calibration curve in different modes
- editing of calibration results (marking of outliers, change of calibration mode)
- creation and print of protocols
  - o calibrations
  - o sample series results
  - o occupancy of the sample tray
- transfer of protocols to pdf- files
- evaluation of control sample results
- archiving of sample and calibration results
- customized presentation of results by means of filter and sorting functions

## General Requirements

hardware: standard PC (min. 4 GB main memory)  
 software: operation system Windows10  
 operation: typical windows user interface

